

# Release Summary

## Medtech Evolution – General Practice

### Version 1.9.0.377 (August 2018)

### Urgent Evolution Update

This release summary provides an overview of the items included within the Medtech Evolution – General Practice Version 1.9.0.377 release:

#### Summary of Resolved Issues

Description of Change	Category
<p>This upgrade fixes an issue where HL7 2.3 RSD formatting characters/tags (example: \H, \N) were not being identified and managed correctly for messages imported into Evolution.</p> <p>This could result in characters and/or words not being displayed fully or correctly within the inbox.</p>	<p><b>Clinical</b></p>

For further information, or any other queries regarding this release, please contact Medtech Support via:

- Medtech Evolution application [**Help ► Contact Support**]
- Online Chat: [www.medtechglobal.com](http://www.medtechglobal.com)
- Email: [support@medtechglobal.com](mailto:support@medtechglobal.com)
- Insight: Register for the Insight Customer Portal and Log a Support Ticket at <https://insight.medtechglobal.com/>