

Release Summary

Medtech Evolution – General Practice

Version 1.9.0.426 (August 2018)

Urgent Evolution Update

This release summary provides an overview of the items included within the Medtech Evolution – General Practice Version 1.9.0.426 release:

Summary of Resolved Issues

| Description of Change | Category |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------|
| <p>This upgrade fixes an issue where HL7 2.3 RSD formatting characters/tags (example: \H, \N) were not being identified and managed correctly for messages imported into Evolution.</p> <p>This could result in characters and/or words not being displayed fully or correctly within the inbox.</p> | <p>Clinical</p> |

For further information, or any other queries regarding this release, please contact Medtech Support via:

- Medtech Evolution application [**Help ► Contact Support**]
- Online Chat: www.medtechglobal.com
- Email: support@medtechglobal.com
- Insight: Register for the Insight Customer Portal and Log a Support Ticket at <https://insight.medtechglobal.com/>